



Search...

RESOURCES (HTTPS://INSIGHT.VISIONSOURCE.COM/RESOURCES)

PURE AND CLEAN (HTTPS://INSIGHT.VISIONSOURCE.COM/PURE-AND-CLEAN)

Pure & Clean



(https://insight.visionsource.com/Vendors/Vendor-Details/vid/538)

What is pure&clean LID PREP? pure&clean LID PREP uses stable electrolyzed Hypochlorous Acid (HOCl). pure&clean LID PREP is an FDA-cleared medical device that is great for irrigating and cleansing the eyelids and lashes while also helping to manage the intact skin around the eye.

pure&clean LID PREP Difference Bottomline: It is safe, powerful, and effective.

Safety First pure&clean LID PREP has shown no side effects and does not contain any surfactants or other irritating ingredients.

Effective, Longer Other HOCl products may only last 30 days after opening, but pure&clean LID PREP is effective for 24 months.

pure&clean LID PREP

- A proprietary stabilized solution of HOCl
- An FDA-cleared 510(K) medical device
- Contains a 0.01% hypochlorous acid formula
- Hypochlorous acid is known to be safe, broad-spectrum and fast acting in solution.
- Hypochlorous Acid is a molecule naturally produced by the human body's immune system.

- Hypochlorous acid cleanses bacteria from around the eyes showing symptoms of blepharitis meibomitis, chronic dry eye, and MGD.*
- Studies have shown HOCl may reduce inflammation, giving patients relief from irritations associated with dry eye

Use and Benefits

- Antibiotic-free and steroid free
- The natural alternative to harsh-chemical products. Essential for cleansing around the eye
- Non-Cytotoxic
- Non-sensitizing and non-irritating to skin
- No sting formula
- No damaging effect on human tissue
- No need to wash after use
- Use two to three times daily for optimum eyelid and eyelash management.

Microorganism Efficacy Test

Test Substance Efficacy Results: These are the results of how the product performed against each test microorganism		
MICROORGANISM	TIME POINT	% REDUCTION
<i>Acinetobacter baumannii</i> ATCC 19606	15 seconds	100%
<i>Enterococcus faecalis</i> (VRE) ATCC 51299	15 seconds	100%
<i>Listeria monocytogens</i> ATCC 15313	60 seconds	>99.96%
<i>Propionibacterium</i> <i>Acnes</i> ATCC 6919	60 seconds	>99.99%
<i>Candida albicans</i> ATCC 10231	60 seconds	>99.9999%
<i>Proteus mirabilis</i> ATCC 7002	15 seconds	>99.999%
<i>Serratia marcescens</i> ATCC 14756)	15 seconds	>99.999%
<i>Staphylococcus aureus</i> (MRSA) ATCC 33592	15 seconds	>99.999%
<i>E. Coli</i> ATCC 8739	60 seconds	>99.999%
<i>Pseudomonas aeruginosa</i> ATCC 9027	60 seconds	>99.999%
<i>S. aureus</i> ATCC 6538	60 seconds	>99.999%
<i>A brasiliensis</i> ATCC 16404	60 seconds	>99.999%
<i>K. pneumonia</i> (CRE) ATCC BAA-2146	60 seconds	>99.9994%
<i>C. Difficile</i> ATCC 43598	60 seconds	>99.9999%
<i>Adenovirus 2</i> ATCC VR-846	60 seconds	>99.99%

All microorganisms were killed in the 15 second contact time. 15 seconds was the quickest time to put neutralizing agents on test cultures.

*Reduction in microbial growth in solution has not been shown to correlate with a reduction of infection in patients. No clinical studies have been performed to evaluate reduction of infection.

For more info click here (<https://insight.visionsource.com/Vendors/ Vendor-Details/vid/538>) to visit the pure&clean vendor page

ORDER NOW ([HTTPS://API.PUREANDCLEANVISION.COM?TOKEN=BEA6A611-2F18-4BF2-BC3F-F3ACCFD5EF1F](https://api.pureandcleanvision.com?TOKEN=BEA6A611-2F18-4BF2-BC3F-F3ACCFD5EF1F))

**Pure and Clean Staff Training
(<http://player.vimeo.com/wmode=opaque>)**

In this webinar, Dr. Pete Kehoe focuses on the benefits and how to effectively implement the Pure & Clean family of products within your office. Feel free to use this video at your next staff meeting to explain the proper use and benefits of these safe, effective and profitable products.



[News + Happenings \(https://insight.visionsource.com/Resources/News-Happenings\)](https://insight.visionsource.com/Resources/News-Happenings)

[New Member \(https://insight.visionsource.com/Resources/New-Member\)](https://insight.visionsource.com/Resources/New-Member)

[Mentor OD \(https://insight.visionsource.com/Resources/Mentor-OD\)](https://insight.visionsource.com/Resources/Mentor-OD)

[Giving Back \(https://insight.visionsource.com/Resources/Giving-Back\)](https://insight.visionsource.com/Resources/Giving-Back)

[Team \(https://insight.visionsource.com/Resources/Team\)](https://insight.visionsource.com/Resources/Team)

[Refer a Friend \(https://insight.visionsource.com/Resources/Refer-a-Friend\)](https://insight.visionsource.com/Resources/Refer-a-Friend)

[Vision Source Next \(https://insight.visionsource.com/VSNext/Home\)](https://insight.visionsource.com/VSNext/Home)

[Pure and Clean \(https://insight.visionsource.com/Pure-and-Clean\)](https://insight.visionsource.com/Pure-and-Clean)

[Member Idea Videos \(https://insight.visionsource.com/Resources/Member-Idea-Videos\)](https://insight.visionsource.com/Resources/Member-Idea-Videos)

© Copyright 2020, Vision Source All rights reserved.

[Terms of Use \(/terms-of-use\)](/terms-of-use) - [Privacy \(/Privacy-statement\)](/Privacy-statement) - [Rules of Engagement \(/Rules-of-Engagement\)](/Rules-of-Engagement)

